

MiRXES BlitzAmp qPCR Master Mix

To Know is to Act

Compatible with most Real-time qPCR machines

With different concoctions (i.e., low/high ROX) available, MiRXES BlitzAmp qPCR Master Mix is designed to support most real-time qPCR machines in the market.

Table 2. Covers the Real-time qPCR instruments that are compatible with BlitzAmp qPCR Master Mix, their ROX requirement and the catalogue numbers of the qPCR Master Mix for ordering.

Real-time PCR Instrument	ROX requirement	BlitzAmp qPCR Master Mix (Rxn)
<ul style="list-style-type: none"> • Bio-Rad® iQ™5, CFX96, CFX384, Opticon • Roche Lightcycler® • Qiagen Rotor-Gene™ • Eppendorf Mastercycler® • Cepheid® SmartCycler® • MIC qPCR 	Not required	1204202 (200) 1204204 (800) 1204205 (1200) 1204206 (2400)
<ul style="list-style-type: none"> • Applied Biosystems® 7500, QuantStudio™, • ViiA7™ • Agilent Mx™ 	Low ROX	
<ul style="list-style-type: none"> • Applied Biosystems® 7000, 7300, 7700, 7900HT, StepOne™, StepOnePlus™ 	High ROX	1204212 (200) 1204214 (800) 1204215 (1200) 1204216 (2400)